

AUTOMATED EXTERNAL DEFIBRILLATOR

DefiMax AED



Intuitive and easy to use semi-automatic defibrillator.

Defibrillator AED dedicated to be used in public access areas, such as shopping malls, sports facilities, workplaces, airports, railways stations, etc., fully compliant with directions of the European Resuscitation Council regarding the use of automated defibrillators.

Thanks to the clear voice messages and graphic instructions, as well as built-in metronome, even a person without a medical training is able to resuscitate a victim with a sudden cardiac arrest until an arrival of paramedics.

Easiness of use, intuitive interface and biphasic, advanced technology ensure the highest reliability of the device.

**EUROPEAN PRODUCT
SOLD AROUND THE WORLD**

Over 25 years of experience, thousands of patient monitors and defibrillators sold in Europe, Asia, Australia, Africa and South America.



EUROPEAN QUALITY



EASY TO USE



RELIABLE

TECHNICAL PARAMETERS

GENERAL PARAMETERS

Power supply	12V battery
Safety class	BF
Safety requirements	MDD: 93/42/EEC, 2007/47/EC; EN 60601-1, EN 60601-1-2, EN 60601-2-4
Weight	3.1 kg without electrodes
Dimensions	275 x 330 x 120 mm

MODES

AED mode

DEFIBRILLATOR MODULE

Electrodes type	single patient use
Impulse type	biphasic with patient's impedance compensation
Energy range	adults up to 300 J children up to 70 J

ADDITIONAL FUNCTIONALITY

CPR metronome
Self-test
Visual battery check

ACCESSORIES

Single patient use electrodes for children
Single patient use electrodes for adults
Battery pack
Transport bag

